





How to

Read AbiBird Activity Graphs and Notifications

The AbiBird sensor monitors movement in a room and compares it to the daily routine set for the sensor in the AbiBird spp. Every hour AbiBird sends the data it has collected to the AbiBird app on your phone. If AbiBird detects an anomaly (movement not matching the set daily routine), it marks the activity as unusual and sends push and email notifications. The sensor also detects vibrations and will send an immediate 'Call for Help' if it detects a strong, continuous vibration.

AbiBird has three states:

- 1. Everything is fine
- 2. Unusual Activity
- 3. Call for Help

Everything is fine!

The *Everything is fine!* state means that the activity in the room matches the routine set in the AbiBird app.

Notifications	Homepage (Sensor List)	Sensor's graph page	
No notifications are sent.	The AbiBird homepage shows an Everything is fine! message for the sensor. My sensors Living Room Everything is fine!	l ''	The example shows there was activity (150 motion events) in the living room between 6:45am and 7:45am. Because the daily routine is set to expect this activity in the morning, AbiBird reports that everything is fine. The cared-for person is up and about as usual.

This guide is for all Smartphone users.



Unusual activity

The *Unusual activity* state means that the activity in the room doesn't match the routine set in the AbiBird app.

Homepage **Notifications** Sensor's graph page (Sensor List) The AbiBird homepage AbiBird marks the A red dot appears on the sensor's graph indicating either; activity as unusual and shows a *Please check* there has been no activity in the room when there should sends a phone and on this person message be; or there's activity in the room when there shouldn't be. email notification. for the sensor. This example shows there was no activity (0 motion events) in the living room My sensors between 7:32am and 8:32am. AbiBird has marked this as unusual Please check in on this person because the daily routine is set for activity support AbiBird notification: Living Room has dete... to begin at 6:30am and alert sensitivity is set for 2 hours. Thus, it's 2 hours since activity was expected to start in the room ■ : and there's been no movement. AbiBird notification: Living Room has detected unusual activity You may want to check on this person as they should be up and about by now. support@stscloud.io NOTE: You can scroll back through the graph to view activity over the last 24 Unusual activity hours. Tap an update (black dot) for the time and number of motion events Unusual activity has been detected by Living Room. You can check the app at any time to see the activity for the day. recorded during that hour. Thank you, This example shows there was activity (24 AbiBird Living Room motion events) in the living room between 10:24am and 11:24am. This message was sent to you by AbiBird AbiBird marked this as unusual because: a) the daily routine is set for activity to end at 10:00am when the person is expected to go out for the day; and b) alert sensitivity is set for 2 hours. Tapping the previous update (black dot) shows there were 3 motion events in the previous hour. Thus, there have been 27 movements over the last 2 hours, which suggests the person stayed home. You may want to call to find out why.



This guide is for all Smartphone users.



Call for Help

The Call for Help! state means that the cared-for person has activated the sensor's vibration sensor.

Homepage **Notifications** Sensor's graph page (Sensor List) AbiBird marks the The AbiBird homepage A yellow dot appears on the sensor's graph indicating a shows a Please check activity as unusual and 'Call for Help' has been triggered. immediately sends on this person message The example shows there was a call for Living Room a phone and email for the sensor. help at 5:26pm. notification. While AbiBird isn't due to send an update My sensors until 5:30pm, a Call for help overrides the Abilitid A call for help has been sent from Living Room hourly update and is marked immediately. Please check in on this person. You would call this person immediately to support Call for help from Living Room check they are okay. Call for help from Living Room Call for help! A call for help has been sent from Living Room. This is activated when the sensor is swiftly tapped a few times or there are multiple loud sounds within a few minutes. You should get in touch with the person and see if they need assistance. AbiBird This message was sent to you by AbiBird